

CORSETS	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)	XXL (6)	XXXL (7)	(8)	(9)
A-1, A-1G, A-1L, A-1A, A1B, A-1C, A-2, B-12, M-1, M-2, M-3, M-4, N-4, N-5, B-1	WAIST	60-70	70-80	80-90	90-100	100-110	110-120	120-130		
B-2, B-3, B-4	WAIST	50-60	60-75	75-90	90-110	110-130				
A-3, A-4, A-5, A-7, A-13, X-17	WAIST	60-65	65-75	75-85	85-95	95-105	105-115	115-125		
A-4A, A-4B, A-4C, A-6, A-8, A-9, A-10, Y-29, N-64	WAIST	70-75	75-80	80-85	85-90	90-100	100-110	110-125		
D-2A, Y-10, Y-52, Y-53	WAIST	60-65	65-70	70-80	80-90	90-100	100-110	110-125		
A-1W, A-2W, N-9, N-9A, N-26, N-26A, N-27, N-30A, BİO-02, N-4A	WAIST		60-75	75-100	100-130					
E-2, E-3, E-8, N-29, N-30, M-5, M-6	WAIST	60-70	70-80	80-90	90-100	100-110				
N-2, N-8, N-17	WAIST		60-75	75-90	90-110	110-130				
MB-2A, M-2A	WAIST									
MB-10	WAIST	12-16-	16-21	21-28	28-36	36-44				
ES-01, ES-02, ES-03, ES-04, ES-05, ES-06, ES-07, ES-08,	WAIST	38	40	42	44	46	48	50	52	54
K-1,	WAIST	38-40	40-42	42-44	44-48	48-52				
K-3, K-3A	WAIST	32-38	38-44	44-50	50-56	56-62				
E-7, E-7A	WAIST	32-34	34-38	38-42	42-46	46-50	50-54	54-58		



KNEE BRACES	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)	XXL (6)	XXXL (7)
V-4A	KNEE		26-32	32-37	37-42	42-46		
R-1, R-2	KNEE		28-30	30-33	33-36	36-39	39-42	
R-3, R-4	KNEE		28-30	30-33	33-36	36-39		
N-31, N-31W, N-32, N-32W, N-33, N-34, N-35, N-35A, N-36, N-37, B-5, X-22, X-23, X-24, X-25, X-45	KNEE	28-30	30-33	33-36	36-39	39-42	42-45	45-48
N-57, N-63, N-63A	KNEE	24-27	27-30	30-34	34-39	39-45		



CERVICAL SUPPORTS	CM.	XXS (0)	XS (1)	S (2)	M (3)	L (4)	XL (5)
C-1, C-5	NECK		22-27	27-32	32-38	38-44	44-50
C-10, C-13, C-13A, C-14, C-15, C-16, C-17	NECK	23-26	26-29	29-32	32-36	36-41	41-46
C-3, C-3A, C-11, C-11B, C-11C	NECK		26-29	29-32	32-36	36-41	41-46
C-11E, C-10E, C-5A	NECK			26-32	32-38	38-46	
C-10R	NECK		22-26	26-32	32-40	40-48	
B-9	NECK		32-36	36-38	38-40	40-43	43-46



WRIST BANDAGES	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)
O-6, B-11, X-5, N-42, N-43, N-44, M-8	WRIST	13-14	14-16	16-18	18-20	20-22
O-11, O-11A, O-11E, O-11C, O-12, N-11, N-43A, N-43B, N-43C, N-43D, N-44A, M-7	WRIST	13-14	14-16	16-18	18-20	20-22
O-11R	WRIST	13-14	14-16	16-18	18-20	
K-11, K-11A, K-11E	WRIST		13-17	17-21	21-25	





ELBOW SUPPORTS	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)		
N-45, N-46, N-47, N-48, N-48A								
B-6, O-7, W-3, R-6	ELBOW	21-23	23-25	25-28	28-31	31-35		
V-5E	ELBOW		21-24	24-28	28-31	31-35		
V-21	ELBOW		22-25	25-27	27-29	29-31		



ANKLE BANDAGES	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)		
B-10, O-10, V-5B, V-5C, V-5D, N-49, N-50, N-51, N-52, N-53, N53A N-53D	ANKLE	18-20	20-22	22-24	24-26	26-28		
V-5A, R-5	ANKLE		18-20	20-23	23-25	25-28		
N-53W, N-53T, D-8E	ANKLE		33-37	37-40	40-45			



RIB BANDAGES	CM.	XXS (0)	XS (1)	S (2)	M (3)	L (4)	XL (5)	
O-25, O-37, O-38	CHEST	55-65	65-75	70-85	80-95	90-105	100-120	
N-1, N-54	CHEST		60-75	75-90	90-105	105-120	RIGHT/ LEFT	
A-12, A-12A, B-8	CHEST			60-80	80-100	100-120		



CALF/THIGH BANDAGES	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)	XXL (6)	
N-38, B-7	THIGH	46-48	48-50	50-52	52-54	54-56	56-60	
N-39	CALF	31-33	33-35	35-37	37-39	39-41	41-43	



VARICOSIS STOCKINGS	CM.	XS (1)	S (2)	M (3)	L (4)	XL (5)		
V-6... V-19 all models	CALF	30-32	32-34	34-36	36-40	40-44		

WASHING INSTRUCTION

Hand wash at 30 C is recommended for all products.

Please close the velcro bands and remove the splints from open pockets before you wash.

Wash slowly in luke-warm water (30 C recommended) frothed with ordinary soap without rubbing. Rinse out in cold water without bending and dehydrate by pressing slowly. Do not dry them under the sun, by the stove or with the iron. Bleaching liquid is not used. Do not iron and do not give them to dry cleaning.