



R-1E KNITTED KNEE SUPPORT

It is knitted from anti-sweating elastic thread. The knee pads are given the anatomical shape of the knee during weaving. It is manufactured in 6 different sizes.

Indication: Used in knee sprains, fiber ruptures, after trauma and surgery, and in sports activities.

R-2E KNITTED PATELLA SUPPORTED KNEE SUPPORT

The knee brace, knitted from anti-sweating elastic rope, is shaped according to the knee anatomy. Extra support is used in the area of the knee brace that comes to the patella. The knee brace is manufactured in 6 different sizes.

Indication: Used in knee traumas, inflammatory events of the knee joint, knee osteoarthritis, fixation of the patella, post-operation and sports activities.



R-3E KNITTED FLEXIBLE UNDERWIRE PATELLA SUPPORTED KNEE SUPPORT

The knee brace, knitted from anti-sweating elastic rope, has been shaped in accordance with the knee anatomy. Extra support has been sewn to the area of the knee brace that comes to the patella, and support has been provided to the lateral and medial ligaments by placing underwires in the columns sewn to the sides of the knee brace. The knee brace is manufactured in 6 different sizes.

Indication: Used in medial and lateral collateral ligament instability and strains, knee traumas, inflammatory events of the knee joint, knee osteoarthritis, fixation of the patella, post-operation and sports activities.





R-4E KNITTED FLEXIBLE UNDERWIND TAPE KNEE SUPPORT

The knee brace, knitted from anti-sweating elastic rope, has been shaped in accordance with the knee anatomy. Extra support has been sewn to the area of the knee brace that comes to the patella, and support has been provided to the lateral and medial ligaments by placing underwires inside the columns sewn to the sides of the knee brace. In order to prevent the knee brace from slipping from the knee, two velcrobands have been placed on the top and bottom to wrap the knee from end to end. The knee brace is manufactured in 6 different sizes.

Indication: Used in medial and lateral collateral ligament instability and strains, knee traumas, inflammatory events of the knee joint, knee osteoarthritis, fixation of the patella, post-operation and sports activities.

R-4EM KNITTED HINGED BANDED KNEE SUPPORT

The knee brace, knitted from anti-sweating elastic rope, has been shaped in accordance with the knee anatomy. Extra support has been sewn to the area of the knee brace that comes to the patella, and metal hinges pressed into the columns sewn to the sides of the knee brace as standard have been placed. In order to prevent the knee brace from slipping from the knee, two velcrobands, one above and one below, have been placed to wrap the knee from end to end. The knee brace is manufactured in 6 different sizes.

Indication: Used in medial and lateral collateral ligament instability and strains, knee traumas, inflammatory events of the knee joint, knee osteoarthritis, knee instability (genu valgus and varum deformity), after arthroscopy.



R-5E KNITTED MALLEOL SUPPORTED ANKLE BRACE

The wristband, which is knitted from anti-sweating elastic rope, provides support to the ankle with special anatomical silicone pads located on the medial and lateral sides, relieves edema, accelerates circulation and provides protection to the ankle. The wristband is produced in 5 different sizes.

Indications: Used in acute and chronic ankle sprains, ligament instability, soft tissue traumas, and during sports activities.



R-5EB KNITTED ANKLE BRACE

It is knitted from anti-sweating elastic rope. During production, the wristbands are given the anatomical shape of the wrist to ensure full support. The wristband is produced in 5 different sizes.

Indications: Used in ankle sprains and ligament ruptures.

R-5EP KNITTED ANKLE BRACE WITH 8 BANDAGES

It is knitted from anti-sweating elastic rope. During production, the wristbands are given the anatomical shape of the wrist to ensure full support. The bandage wrapped around the wrist in the shape of 8 provides extra support. The wristband is produced in 5 different sizes.

Indications: Used in ankle sprains and ligament ruptures.



R-6E KNITTED SILICONE SUPPORTED ELBOW BRACE

The elbow pad, knitted from anti-sweating elastic thread, is designed in accordance with the elbow anatomy. It provides extra compression with silicone supports sewn to the area where lateral and medial epicondylitis occurs. It is produced in 5 different sizes.

Indications: Used in elbow traumas, elbow lateral and medial epicondylitis, bursitis, elbow osteoarthritis, ligament injuries, post-operative recovery process and after arthroscopy.

R-6EB KNITTED ELBOW BRACE

The elbow pad, knitted from anti-sweating elastic thread, is designed in accordance with the elbow anatomy. It is produced in 5 different sizes.

Indications: Used as a protector in elbow traumas, elbow osteoarthritis and arthritis, ligament injuries, post-operative recovery process and sports activities.



R-7E KNITTED LUMBOSACRAL CORSET

The back of the corset, which is knitted from anti-sweating elastic rope, is supported by 3 steel underwires and the corset is manufactured in 5 different sizes to suit waist measurements. The corset, which has finger holes at the front so that the patient can put on the corset by himself, is produced in gray.

Indications: Used in acute and chronic lumbago, lumbosacral traumas, lumbar hernias, lumbar spondylosis, spondylolisthesis, post-operative conditions.



R-8E KNITTED WRIST SPLINT

The wristband, woven from anti-sweating elastic rope, is designed to suit the wrist anatomy. Full stabilization is provided thanks to the aluminum underwire located in the palm. It is supported by a velcro from the top. It is compatible with the right or left wrist. It is produced in 5 different sizes.

Indications: Used to keep the hand in a resting position in carpal tunnel syndrome, hand-wrist traumas, sprains, injuries and after surgery.





R-9E KNITTED KNEE BRACE

It is knitted from anti-sweating phosphorescent elastic thread. The knee pads are given the anatomical shape of the knee during weaving. It is manufactured in 4 different sizes.

Indication: Used in knee sprains, fiber ruptures, after trauma and surgery, and in sports activities.

R-10E KNITTED PHOSPHOROUS ANKLE BRACE

It is knitted from anti-sweating phosphorescent elastic thread. During production, the wristbands are given the anatomical shape of the wrist to ensure full support. The wristband is produced in 4 different sizes.

Indications: Used in ankle sprains and ligament ruptures.



R-11E KNITTED PHOSPHOROUS ELBOW BRACE

The elbow pad, knitted from anti-sweating phosphorescent elastic thread, is designed in accordance with the elbow anatomy. It is produced in 4 different sizes.

Indications: Used as a protector in elbow traumas, elbow osteoarthritis and arthritis, ligament injuries, post-operative recovery process and sports activities.

R-12E KNITTED PHOSPHOROUS WRISTBAND

The wristband, which is knitted from anti-sweating elastic phosphorescent thread, is designed to suit the wrist anatomy. It is produced in 4 different sizes.

Indications: Used as a protector in hand-wrist traumas, sprains, strains and wrist pains, and in sports activities.



R-12P KNITTED WRIST BAND WITH PHOSPHOROUS BANDAGE

The wristband, which is knitted from anti-sweating elastic phosphorescent thread, is designed to suit the wrist anatomy. The bandage on the wrist provides extra support. It is produced in 4 different sizes.

Indications: Used as a protector in hand-wrist traumas, sprains, strains and wrist pains, and in sports activities.